## Added features in Release RSVPKeyboard2013.4.00:

- Images in the presentation, will be shown in the same color as characters (according to the setting in the parameter files).
- It is possible to save "daq" files in "mat" format and use them in "noAmps" mode or "offlineAnalysis" to do this in the command window please type the following commands:
  - addpath(genpath('.'))
  - 2. loadSessionData

Then you can choose the appropriate "daq" file and it will be saved in "mat" format at the same location with the same name. Finally, instead of loading the "daq" file in any operation you can choose the corresponding "mat" file.

- In the "userParameters" in simulation section there are three new options :
  - 1. RSVPKeyboardParams.Simulation.ExportToExcel.Enabled = 1;
  - RSVPKeyboardParams.Simulation.ExportToExcel.Filename = [];
  - RSVPKeyboardParams.Simulation.ExportToExcel.Folder = [];

Using these options you can choose to have the simulation results exported to and excel file. "FileName" and "Folder" can be left empty which in that case file will be saved using the default file name and at the corresponding session directory.